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Hospital Universitario de Niguarda, Milan, Italy. October 2022

Cardiologic Institute of Galle. Sri Lanka. 2014

Pediatric Hospital of Ho Chi Minh City, Vietnam. June 2023

Surgical Hospital of Ashgabat, Turkmenistán, April 2018

Civil Hospital of Guadalajara – Mexico. November 2018

Indira Gandhi Memorial Hospital. Male – Maldives. October 2023

Provincial Hospital of Bengo, Angola. February 2024

National Nurse Congres. Nanning, China. August 2019

Military Hospital of Alexandria. Egypt. September 2015

CEDIMAT. Dominican Republic 2016


NSW Congress – Coffs Harbour. Australia. March 2024

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Porquê?

https://www.rki.de/DE/Content/infekt/Krankenhaushygiene/Kommission/Downloads/Hygiene_Requirements_Medical_Devices_2012.pdf?__blob=publicationFile
Table 1. Risk Assessment and Classification of Medical Devices before Reprocessing

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Overcoming the problem of residual microbial contamination in dental suction units left by conventional disinfection using novel single component suction handpieces in combination with automated flood disinfection .Boyle MA, et al J Dent. 2015 Oct;43(10):1268-79

“A proposed cleaning classification system for reusable medical devices to complement the Spaulding classification”

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Kremer T, Rowan NJ, McDonnell G, A proposed cleaning classification system for reusable medical devices to complement the Spaulding classification, *Journal of Hospital Infection*, <https://doi.org/10.1016/j.jhin.2023.11.018>. 22 November 2023

CONGRESSO NACIONAL^{EM} **REPROCESSAMENTO** DE **DISPOSITIVOS MÉDICOS**
SANTO TIROSO | PORTUGAL

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ISO 17664-1:2021(E)

Annex C
(informative)

Classification of medical devices

C.3 Medical device design groups for processing

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Canulado Sujo
<https://www.inspektor.com/videos-imagens/>

Canal endoscópico danificado
<https://www.inspektor.com/videos-imagens/>

Interior Lumen AFTER Sterilization

Interior Lumen AFTER Sterilization

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CyberKnife® S7™


« THE EVOLUTION AND FUTURE OF »
ROBOT-ASSISTED SURGERY

THE ROBOTIC SURGERY MARKET IS EXPECTED TO SURPASS
\$7 BILLION BY 2025.


<https://online.regiscollege.edu/blog/robot-assisted-surgery/>
<https://mantadesign.com/the-future-of-robotic-surgery/>

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
Novel concept of cleanliness of instruments for robotic surgery.

<http://dx.doi.org/10.1016/j.jhin.2016.04.009>

Journal of Hospital Infection 93 (2016) 360–361
Available online at www.elsevier.com/locate/jhin
Journal of Hospital Infection
Journal homepage: www.elsevier.com/locate/jhin

Opinion
Novel concept of cleanliness of instruments for robotic surgery

Introduction
Robotic systems for surgery have been developed over the last 30 years, and currently use of the da Vinci® Surgical System (Intuitive Surgical Inc., Sunnyvale, CA, USA) predominates in many countries. Robotic systems provide surgeons with many benefits such as refined manipulation of surgical instruments. On the other hand, they bring about difficulty for personnel in central sterile supply departments, because these instruments may be difficult or impossible to clean adequately. This article discusses concerns about the cleanliness of robotic instruments.

Cleaning of robotic instruments
Robotic instruments mainly consist of tips, shafts and housings. They are mounted on robotic arms and manipulated by remote control from a console. Most robotic instruments are semi-reusable and need to be decontaminated after each use up to 10 times. Unlike other endoscopic instruments, robotic instruments are cleaned without disassembly. Even worse, their complex structures also interfere with cleaning and reprocessing.

According to the manufacturers' instructions, cleaning of robotic instruments should be performed by either a meticulous manual procedure or a manual plus mechanical procedure. The manual cleaning procedure includes priming with enzymatic detergent, brushing the outer surface, flushing ports, ultrasonication, and rinsing. As manual brushing is only applied to the outer surface, flushing of shafts and housings is relied upon to remove contamination of the inner surface. A major concern about manufacturers' instructions is that they do not provide either benchmarks or evaluation methods for cleaning efficacy. Manufacturers probably have not established standards for the efficacy of mechanical cleaning either, although several washers have already been approved.

Residual contamination after cleaning
In our own research on the cleanliness of 24 robotic instruments, they still contained 21.7 mg protein/instrument on average (unpublished data) following the cleaning procedure recommended by the manufacturer. The amount of residual protein was much greater than that on minimally invasive instruments.

What sort of instruments are difficult to clean?
The reason why manufacturers' instructions do not propose relevant methodologies for evaluating cleanliness would partly be a lack of the concept of cleanliness in their design. The structure of robotic instruments is only designed to achieve fine manipulation and durability, whereas cleaning adequacy would not be the top priority in assembling the instruments. There is an opinion that all robotic instruments should be single-use devices in order to prevent cross-contamination via surgical instruments if they cannot be cleaned adequately. However, this so-called 'single-use policy' would greatly increase the cost. Such cost could be reduced if there were methods to determine accurately the cleanliness of instruments. In our view, measurement of protein residue on instruments would be a substitute of the conventional classification based on geometric features. The pattern of protein residue in repeated extraction might be the structure of instruments in our preliminary study, which a method could be utilized to estimate how difficult it is to clean an instrument.

Conclusions
A standard classification should be established to evaluate the cleanliness of robotic instruments based on how difficult it is to clean an instrument. This would also reduce the cost of robotic surgery by avoiding extensive use of single-use devices.

Conflict of interest statement
None declared.

<http://dx.doi.org/10.1016/j.jhin.2016.04.009>
0195-6701/© 2016 The Healthcare Infection Society. Published by Elsevier Ltd. All rights reserved.

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EN ISO 15883-5

Washer-disinfectors - Part 5: Performance requirements and test method criteria for demonstrating cleaning efficacy (ISO 15883-5:2021)

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